

CATALOG

CHOUKROUN PRF™ SYSTEM



PRF PROCESS™
CHOUKROUN

49, rue Gioffredo • 06000 NICE • FRANCE
www.a-prf.com • info@a-prf.com
00.33.493.855.890

THE EQUIPMENT

▶ DUO Quattro Centrifuge:





MD Class IIa - CE 2797

The new centrifuge DUO Quattro has a «push» button that affords you to select between 6 different PRF protocols. The values of these settings are preset but not locked (you can set other values if needed) .

Additional position : Custom program

Manual : Free settings

Simply push the button to select the desired mode of operation, the corresponding LED will light on. The machine is ready. Simply press the button  to start the centrifugation.



6 PROTOCOLS

▶ PRF BoX :

« PRF Box » allow you to get the fibrin membranes rich in platelets always hydrated and of constant thickness, but also to recover the exsudate rich in integrins: Vitronectin and Fibronectin. You can also produce “plugs” of PRF in white teflon cylinders.

- ① Multibacs in teflon, with compartments of different sizes to easily create large membranes or sticky bone.
- ② Stainless steel crusher, to crush PRF clots

BOX



 **New**



PomPac™



The most important innovation in the PRF field in the last 20 years ! In order to constantly improve the technique and results of PRF, discover our new product : PomPac™

PomPac™ will guarantee that our users will always get larger PRF clots and liquid PRF with 50% more volume!



PomCol™



For sampling of 8 or more tubes ...use PomCol™


Procedure: Cool the sampler tubing to prevent blood coagulation in the sampler.



PomSwing



The PomSwing is composed of 1 jar and a Roll. It allows you to shape the Sticky Bone quickly and easily.



PomPress



Anti-stress balls, should be used during the blood drawing in order to relax the patient. Squeezing the ball will also increase the blood flow in the veins.





The instruments :



PRF Scissors

Allows to separate easily the fibrin clot from the red cells and to cut membranes into small pieces.



PRF Forceps

Is used to remove and manipulate membranes.



PRF Forceps giraffe

The forceps giraffe is intended to tunneling the PRF membranes for surgical soft tissues procedures: introduce the forceps in the tunnel, open the forceps, block the PRF membranes and pull to place the membranes in the correct position.



PRF Tray

- Used to prepare membranes and too:
- Separate the fibrin clot from the red blood cells
 - Cut the membranes
 - Fold the membranes
 - Make «punches»




PRF Bowl

Is used to mix the membranes cut into small pieces with biomaterials. It can also be used for the cutting: Put the membrane into the bowl and cut it with PRF scissors.



The instruments :

1
In the Sinus: this plate is used to retrieve the PRF membrane (always) folded twice and slide it into the sinus to the desired location easily. Introduce the pointed end of the pad 1 in the sinus window and drag the PRF membrane on the Schneiderian membrane. Same process for a collagen membrane.



Mini Tray



2

Is used to position the PRF membranes easily, anywhere, using the pad as a cake server. Drag the PRF membrane on Pad 2 and place it on the desired surgical area. Same process for a collagen membrane...

Allows to compact and tamp the bone filling materials in the recipient sites (alveoli, sinus, goshawks implants).

Allows to compact and tamp the bone filling materials in the recipient sites (alveoli, sinus, goshawks implants).

Allow to create and mold easily Sticky Bone

▶ The instruments :

Tube-holder

Is used to hold the tubes and let them rest after the blood drawing



Tourniquet

It must be tightened to the patient's forearm 10 cm above the site of venipuncture.



Blood Collector




Sticker

Is used to maintain in place the collector during the blood drawing



▶ Tubes /



A-PRF tubes

S-PRF tubes

i-PRF tubes



Clot
A-PRF Liquid



PRF Liquid
Sticky Bone
Large
membrane





injectable PRF
Aesthetics
Orthopedics

A-PRF tubes

S-PRF tubes

i-PRF+ tubes





MEMBRANE, CLOT, PLUG

A-PRF 10mL tubes
Medical device Class IIa.

Box of 100 tubes
CE 0425

Made in EUROPE

Dr. Choukroun Glass A-PRF™ Tubes, No Additives (Sterile) 100/Box

- Dr. Choukroun A-PRF™ tubes are used for the production of A-PRF membranes
- Dr. Choukroun A-PRF™ tubes are made of glass with no additives
- Dr. Choukroun A-PRF™ tubes are 10ml
- Dr. Choukroun A-PRF™ tubes are non-toxic & non-pyrogenic
- Increased vacuum ensures that the blood draw is completed quickly
- Packaged 100/box
- Dr. Choukroun A-PRF™ tubes are sterile and individually packaged



STICKY BONE


S-PRF 10mL tubes
Medical device Class IIa.

Box of 24 tubes
CE 0425

Made in EUROPE

Dr. Choukroun Plastic S-PRF™ Tubes

- Dr. Choukroun S-PRF™ tubes are used for the production of plasma for sticky bone
- Dr. Choukroun S-PRF™ tubes are made of plastic with no additives
- Dr. Choukroun S-PRF™ tubes are 10ml
- Dr. Choukroun S-PRF™ tubes are non-toxic & non-pyrogenic
- Increased vacuum ensures that the blood draw is completed quickly
- Micro-texture helps prevent clot formation
- Packaged 24/box
- Dr. Choukroun S-PRF™ tubes are sterile and individually packaged



FABRIQUÉ
EN EURPE
Produced in EUROPE

RESERVED FOR USE :
Aesthetics and Orthopedics






i-PRF tubes 13mL
Medical device 25 tubes

CE 0425



Made in EUROPE

Dr. Choukroun Plastic I-PRF+ Tubes For Injections and Aesthetics

- Dr. Choukroun I-PRF+™ tubes are used for the production of plasma for facial aesthetics & orthopedics
- Dr. Choukroun I-PRF+™ tubes are made of plastic with no additives
- Dr. Choukroun I-PRF+™ tubes are 13ml
- Dr. Choukroun I-PRF™ tubes are non-toxic & non-pyrogenic
- Increased vacuum ensures that the blood draw is completed quickly
- Micro-texture prevents clot formation
- Packaged 24/box
- Dr. Choukroun I-PRF+™ tubes are sterile and individually packaged
- [PRF Tube Quality Comparison](#)

3 tubes are now available :		A-PRF 12 / DUO / Quattro 	PC 02 	
A-PRF™ Red tubes 	Clots Membranes Plugs	1300rpm - 14min	1500rpm - 14min	
S-PRF™ Green tubes 	Sticky bone Big membranes	1300rpm - 14min	1500rpm - 14min	
	PRF Liquid Injectable <small>(S-PRF green tubes replace i-PRF orange tubes)</small>	i-PRF (women)	700 rpm - 3 min	800 rpm - 3 min
i-PRF™+ Purple tubes 	ONLY for use in aesthetics and orthopedics	i-PRF^M (men)	700 rpm - 4 min	800 rpm - 4 min
		i-PRF⁺ (women)	700 rpm - 5 min	800 rpm - 5 min
		i-PRF⁺M (men)	700 rpm - 6 min	800 rpm - 6 min

SAFETY Blood Collection/Infusion Set
Instructions for Use



Production location:
Nipro (Thailand) Corporation Limited,
10/2 Moo 8 Bangnomko, Sena,
Phra Nakhon Si Ayutthaya 13110, Thailand
Made in Thailand
Distributed by Greiner Bio-One GmbH, Austria



Greiner Bio-One GmbH
Bad Haller Straße 32
4550 Kremsmünster
AUSTRIA

www.gbo.com/preanalytics
faq.pa@gbo.com
(+43) 07583/6791



Please visit our Download Center
at www.gbo.com to download
the full instructions for use.
see IFU 980224

POLYSTERIBOX®



The polyesteribox® is a reusable container for sterilization, transport and storage of your instruments. This box is made of «DURADEX» (PPSU), This material allows higher stability, durability and resistance. (up to 150°C)

Filter duration **5 years** or **1250 sterilizations**
Duration of sterilization: **6 months**

Cleaning:

Clean with an alkaline detergent and rinse with water.

Sterilization:

- ✓ Autoclave sterilization up to 134°C.
- ✓ Chemical disinfection at a maximum temperature of 60°C
- ✓ Thermal disinfection at a maximum temperature of 93°C

⚠ Do not use drying aids! Risk of damage! The material becomes opaque.

The polyesteribox® L can contain all the PRF kit for an optimal and practical sterilization.



Caution: Remove the crusher from the PRF BoX™ before sterilization




Arrange all Instruments in the Polyesteribox®



Close, place the sterilization label. Then proceed with the sterilization. After that seal with the safety seals

Delivered with :



Silicone positioning mat



Sterilization seals



Safety seals

Optional :



Strainer



49, rue Gioffredo • 06000 NICE • FRANCE
www.a-prf.com • info@a-prf.com
00.33.493.855.890